

BACtrack integrates a smartphone app and police-grade breathalyzer to bring blood-alcohol content (BAC) wirelessly to smartphones. Features Bluetooth connectivity, innovative tracking, video verification, and violation alerts.

## Technical Specifications

## **BACtrack Specifications**

- *Size:* 1.75" x 2.75" x 0.63"
- Weight: 1.75 oz
- Battery Operation: Up to five (5) days
- Battery Charging: Via standard AC outlet with provided micro USB cable

- Sensor Technology: Xtend™ electro-chemical fuel cell
- BAC Range: 0.000 0.400%
- Required Blow Time: Five (5) seconds
- Warm-up Time: Approx. five (5) seconds
- Bluetooth Connection Time: 60 seconds



## **BACtrack Smartphone App Specifications**

- Compatible Smartphones: Android and iPhone (compatible models vary)
- Application Download: Google Play Store (Android), iTunes Store (iPhone)
- · Login Maintained: After initial enrollment/login, participants cannot log out of the app
- Participant Profile Photo: Photo taken via the app is required for participants to complete enrollment and take tests
- Participant Test Notification: Audio alert is generated by the app to notify the participant to take a test
- Step-by-Step Instructions: Participants are guided through the enrollment and testing processes
- Participant Video Verification: Tests require the participant to video themselves while taking the breath test
- Test Schedules: Random, scheduled or manual tests can be set as a one-off or repeated weekly

